


Lake



River



Qiagen Blood & Tissue, w/ Beads



Qiagen Blood & Tissue, w/o Beads



Break cartridge, Qiagen PowerSoil



Qiagen Blood & Tissue, w/ Beads




Qiagen Blood & Tissue, w/o Beads



Break cartridge, Qiagen PowerSoil

Sterivex™



bcl2fastq



Demultiplex
(Claident)



Amplicon
Sequence Variant
(DADA2)



Taxa assignment
(Claident)



Statistical analysis
(Phyloseq)

Sampling + Filtration

DNA extraction, Library preparation

Sequencing

Data analysis